

NOTEBOOK NO. 46
ISSUED TO Yaxson D. Estrada - Amnarij
ON 4/17/07 2007
DEPARTMENT Media Prep
RETURNED 4/10 2007

SCIENTIFIC NOTEBOOK COMPANY
2831 LAWRENCE AVENUE
STEVENSVILLE, MICHIGAN 49127
(800) 537-3028 - <http://www.snco.com>

TLE Media Prep Training

Project No. _____
Book No. _____

Page No. _____

Objective : To learn the overall activities ^{y.d.} ground of media preparation

Procedures

1. Media Lab Prep PPE

- always wear safety glasses, gloves and lab coat during media preparation.
- dust mask must be worn during preparing media
- Take off the lab coat when leave the room.
- Clean safety glasses by running through water and soap after use.
- Wear close-toe shoes.
- Take extensive awareness when dealing with chemicals.

2. Media ordering process

(1) The requesters must submitted media request form 10 working days before they need.

- Fill in the date requested, the name of requester, date needed (10 working days)
- project or P.I name and account number of the project.
- name of media or reagent
- media recipe book SOP number (RCXXX)
- volume, quantity and other specificity of the media (The media that need to be prepared fresh will be prepared in the morning of the day needed)

- PI signature.

(2) Stamp the form with the date and time and leave the form in the "IN" drop box

(3) Picking-up order

- After the media are prepared, fill in the bottom part of the media request form
- The batch number X.YYYY-Z ; X is the autoclave number.
y is the autoclave ^{year} cycle.
Z is the letter of specify each batch
- After ^{y.d.} sign the requester sign the name and date when pick up the media, the form is put in the "OUT" box.

Assessed & Understood by me,

Amor Blosdyng

Date
1/23/2007

Invented by:
Juwana Lopez

Recorded by:

Date
1/30/07

To Page No. _____

3. Lab rules.

- Let the media lab know you will be working in the lab
- Wear goggles and clean them with soap and put them into place.
- Wash your glasses after use
- Prepare 70% ethanol from 100% ethanol under the hood
- Contact media prep lab technician if there is any problem (when the button at the tank turns red)
- Do not accumulate waste. Put the trash container outside to be picked up.
- Always keep your area clean. Wipe the bench, balance or ^{ID} media spill with 70% ethanol and paper towel.
- Put the equipment into places.
- Clean up anything dropped or spilled on the floor except chemicals, notify technician.
- Ask when have question
- Let the media technician know when supplies are running low.
- Wash your hands before leaving the lab.

4. Walkin Cooler Policy.

- No inoculated plate store in a cooler (any inoculated items not allowed)
- Tubes must be in a secondary container tray, and plates must be in secondary container box.

5. Media Preparation.

- Follow the instruction in the SOP
- media sterilization is usually 121°C 45 minutes, unless other conditions are specify
- measure the ^{pH of} media prepared before and after sterilization.
- In broth, after autoclaving, let it cool on a bench. Also, cover the test tube broth with plastic wrap to prevent evaporation.
- In agar, after autoclaving, put it in a water bath that has been set at 45-47°C. to let cool to 45°C before pour plate.
- After pour plate, let it cool overnight

Witnessed & Understood by me,
Yamun Debata

Date
1/23/2007

Invented by: *Urvana Sena*
Recorded by:

Date
1/30/07

n Page No. _____

Objective: To prepare media agar

Procedures.

1. ~~Look up~~^{Y.D.} ~~locate the~~^{Y.D.} Take out the recipe binder and ~~look up~~^{Y.D.} locate the ingredients.
2. Prepare the ingredients on the table
3. Use the appropriate glassware for the appropriate volume (The media should be less than $\frac{2}{3}$ of the maximum vol. of container)
4. Use cylinder to measure DI water that is needed in preparing the media. Then pour into the glassware prepared. Mark the level of that volume.
5. Discard some of the water
6. Measure the media ingredients needed according to the recipe using weight boats.
7. Carefully pour the measured ingredients into the water.
8. ~~Measure~~^{Y.D.} Use cylinder to fill in the DI water into the glassware to the volume marked.
9. ~~Pull up~~^{Y.D.} Drop a magnetic bar into the flask (the glassware)
10. Pull out some foil and cover the flask.
11. Pull out some autoclave labelling tape and put it on top of the flask. Label the batch letter.
12. Stir and heat mix the agar media, watch it carefully to prevent boiling.
13. Fill in the Daily media batch sheet. (batch sheet binder)
14. Measure the pH of the agar media prepared.
15. Autoclave 121°C , 45 mins, liquid
16. Take out and put it in a waterbath which is set to 45°C .
17. Wait until the media cool down. Then, pour plate. (stir before plating)

Activities

1. Prepare tryptic soy agar 700ml by following the procedures.
2. Prepare sorbitol Macconkey agar by following the procedures.
3. ~~Prepare slant~~^{Y.D.} pipette slants (batch number 3.400 H)

detailed description of how each media was prepared must be included

To Page No. _____

nessed & Understood by me,

Yours *Dilindang*

Date

1/24/2007

Invented by:

Recorded by:

Dilindang Jaya

Date

1/30/07

Project No. _____

Book No. _____

TITLE Tryptic Soy AGAR PREPARATION

From Page No. _____

Objective: To prepare tryptic soy agar, (700 ml)

Procedures

1. ~~Bring out the appropriate glassware and use cylinder to measure.~~ P.D.
1. Bring out the appropriate glassware (flask 1L)
2. Use cylinder to measure 700 ml ~~and~~ and pour into the flask. Then, mark the volume ~~of the glass~~.
3. Discard some or about half of the water
4. Look at the recipe and calculate the amount of agar needed

1 L	suspend	40 g	of	tryptic soy agar.
700 ml	suspend	$\frac{40}{1000} \times 700 = 28$	g	of tryptic soy agar.
5. P.D.
5. Weigh 28 g of tryptic soy agar and pour in the flask
6. Add DI water to the volume needed (700 ml), indicated by the line marked.
7. Drop a magnetic bar into a flask.
8. Cover the flask with foil P.D. and label with autoclave tape.
9. Stir and heat mix the agar media. Watch carefully to prevent boiling
10. Fill in the daily media batch sheet.
11. Pour some media into a plastic bucket to measure pH later.
12. P.D. Sterilize by autoclaving for 121°C, 45 minutes liquid.
13. P.D. Let the media cool down in the 45°C water bath
14. Pour plate using aseptic technique. Stir first
15. Take small sample to test pH after autoclaving

Witnessed & Understood by me,

Yaron Dilschey

Date

8/3/07

Invented by: -

Viriana Lopez

Recorded by:

To Page No. _____

Date

2/8/07

TITLE Sorbitol MacConkey agar preparation

Project No. _____
Book No. _____

From Page No. _____

Objective : To be able to prepare sorbitol MacConkey agar by the standard pouring procedures.

Procedures

1. Prepare 700 ml of MacConkey agar.
2. Bring 1 L flask and use DI water and cylinder to measure the 700 ml scale of the flask
3. Discard some water
4. Look at the recipe and calculate the amount of agar needed.

1L	Suspend	50 g of sorbitol MacConkey agar.
700 ml	Suspend	$\frac{50}{1000} \times 700 = 35 \text{ g}$

5. Weigh 35 g of sorbitol MacConkey agar and carefully pour into the flask.
6. Follow the procedures 6 - 14 of the tryptic soy agar preparation in page 4

Prepared & Understood by me,
Amal Dabheri

Date
2/3/07

Invented by:
Vincent Daza
Recorded by:

Date
2/8/07

To Page No. _____

Project No. _____

Book No. _____

TITLE UVM broth preparation, weekly and monthly housekeep

From Page No. _____

Objective :- To be able to prepare uvm broth
- prepare 5 racks of uvm broth ; each contains 9.5 ml

2) UVM PREPARATION

Procedures.

- 1) Calculate the amount of media needed to prepare.
($5 \times 9.5 \times 9.5$) + 300 = 9.7 L. (make extra 3 tubes for each rack for ^{QC} check)
- 2) Bring out the appropriate glassware and use ^{DI} cylinder to measure DI water and mark the right volume.
- 3) Discard some of the water
- 4) Weigh the correct amount of agar needed (52 g uvm media per 1 L water)
- 5) carefully pour into the glass container
- 6) Use cylinder to fill in DI water to the volume marked.
- 7) Use magnetic bar to help mixing the media under heat.
- 8) After mixing the media thoroughly, transfer into tubes by using Omnispens ^{Dis} - Dispensers.

8.1) Operating Omnisense Dispenser.

- a) wearing clean gloves turn on the pump (from the back)
- b) manual should come up on the view screen, press auto/MAN until "Automatic" is shown.
- c) press SEL 9 times until Recall ~~X~~ comes up and press the program you want (Program 3: 9.5 ml). Then, press ENTER.
On the screen you should now see the volume you selected
If you press the wrong program press SEL twice until you see the "RECALL" on the screen. Press CLEAR once followed by the program you want then ENTER
- d) If you ^{did} pressed the right program and you can verify it by looking at the volume shown on the screen press AUTO/MAN 3 times until you see AUTOMATIC on the screen.
- e) Press START when you are ready to begin and STOP when you are done

8.2) Prepare 2 flasks :- one for DI water and one for ^{dis} waste disposal

8.3) The tube with the glass tip will take the water in and the tube with plastic tip will be the route that water comes out.
(This procedure is to clean the dispenser's tubes) \approx 150 ml of water

8.4) Lift the water-in tube to let all the water in the tube comes out

8.5) Dip the ^{water} tube with the glass tip in the media broth. Click start and let it run once into the ^{media} media container before dip it into the glass tubes, keep some media for pH measurement.

Witnessed & Understood by me,

Jan Dalang

Date
2/4/09

Invented by:

Recorded by:

Date

To Page No. _____

Project No. _____

Book No. _____

Page No. _____

- 2.6) Cap the glass tubes. Not pushed down.
- 2.7) After finish, run with water again (2000 ml)
- 2.8) Turn off the pump.
- 2.9) Clean the bench and area with EtOH.
- 3) Label Y.D. Put the auto-clave tape on top of the rack and label with letters.
- 4) Sterilize by autoclave at 121°C 45 mins
- 5) Cool on the bench overnight. Label the batch number on each rack. Push down caps
- 6) Wrap the racks with plastic wrap.
- 7) Label Y.D. Put a labelling tape on the top.
- 8) Store it in a cooler.
- 9) Weekly housekeeping.
- 1) Clean out side the water bath with EtOH
- 2) Clean the pH meter and balances with EtOH
- 10) Monthly housekeeping
- 1) Take out all the water in water bath. Clean with EtOH. Fill in clean DI water.
- 2) Dust the shelves and cabinets, desks.
- 3) Check the media in the refrigerator for expired media.

To Page No. _____

Prepared & Understood by me,

Rishu Dhillon

Date

2/4/07

Invented by: -

Virendra Singh

Recorded by:

Date

2/8/07

Project No. _____

Book No. _____

TITLE 0.01% peptone broth preparation.

From Page No. _____

Objective: Prepare 4 l of 0.01% peptone broth (2.5 ml tubes) ~

Procedures

- 1) Calculate the amount of media needed to prepare
(5 racks \times 75 \times 9.5 ml) + 300 = 3.8 L. \approx 4 L.
- 2) Bring out the appropriate glassware and marked the right volume by using cylinder and DI water to measure the water.
- 3) Discard some water
- 4) Weigh 2.4 g into 4 l of water.
- 5) Follow the UVM procedures from no. 5 ^{exp} see page 6)

To Page N

Witnessed & Understood by me,

gaurav Patel

Date

2/4/07

Invented by:

Urvashi Singh

Recorded by:

Date


2/8/07

Objective : Prepare 5x700 ml of BHI A ^{Y.D} plates.

Procedures

- 1) Bring of 2 L flasks and marked the point ^{Y.D} where 700 ml lies.
- 2) Discard some amount of water.
- 3) Calculate the amount of agar needed. ($\frac{1000 \text{ ml}}{700 \text{ ml}} \rightarrow 52 \text{ g}$)
- 4) weight 36.4 g of BHI A ^{Y.D} using neighbours and carefully pour it into the water
- 5) Fill in the water to ^{Y.D} the 700 ml in volume (to the marked line)
- 6) Drop in the magnetic bar, cover with foil, label the name of agar and batch letter on the autoclave tape and put it on the ~~the~~ covered foil.
- 7) Heat mix until the agar melt and homogenize, pour some agar to measure pH.
- 8) Autoclave 121°C, 45 mins, liquid.
- 9) ~~put in~~ ^{Y.D} Cool down in 50°C waterbath. (around 1 hr.)
- 10) pour plate using aseptic technique, and measure the pH of ^{Y.D} the each batch.
- 11) ~~Cool down~~ ^{Y.D} Let agar solidify overnight

Clean the pH meter, waterbath, and balances
 Clean Desks, bench, laminar flow, waterbath in room 314.



Checked & Understood by me,
Amir Hilman

Date
2/6/07

Invented by:
Viviana Lopez

Date
2/8/07

Recorded by:

Project No. _____

Book No. _____

TITLE ^{Y.P.} ~~Buffer peptone water preparation~~ BHA preparation

From Page No. _____

~~Objective: To prepare 3 racks of 1 ml of buffer peptone water broth.~~

~~Procedures~~

Objective: To prepare 5 x 100 ml Brain Heart Infusion Agar.

Procedures.

- 1) Marked the glasswares used at the level 100 ml, using DI water.
- 2) Pour out some water ^{Y.P.}
- 3) Weight 36.4 g of BHA (100.715) in each 100 ml (since 1 L need 32 g)
- 4) Add water to 100 ml as marked on the glasswares.
- 5) Cover ^{Y.P.}. Drop a stir bar in, and cover with aluminium foil, labeling tape. Then, Heat mix,
- 6) Pour some media to take pH later.
- 7) Autoclave liquid 45 min, 121°C
- 8) Put in a waterbath (48°C) Wait for cool down 40 mins.
- 9) Pour plate - and separate some of the media to measure pH later.
- 10) leave the ~~agar~~ agar to cool down.

*Note: Also include the other activities you do: Kill, housekeeping, etc To Page No. _____

Witnessed & Understood by me, *gauri Deshpande*

Date 2/9/07

Invented by: *Urvana Jena*

Date 2/15/07

Recorded by:

TITLE BHIA and peptone water preparation.

Project No. _____

Book No. _____

From Page No. _____

Objective: To prepare 5 x 700 ml BHIA and 3 racks of 9 ml.

I.D. (I) Brain Heart Infusion Agar. 5 x 700 ml.

Procedure.

- 1) Use cylinder to measure 700 ml DI water and mark the flask.
- 2) Pour some DI water out. ^{I.D.}
- 3) Weight 36.4 g of media in to each ~~700 ml~~ flask.
- 4) Pour ~~some~~ ^{I.D.} DI water until reach 700 ml.
- 5) Cover with foil and heat mix.
- 6) Pour some media out to measure pH before autoclave.
- 7) Autoclave 121°C 45 min, ^{I.D.} lig.
- 8) After taking out from the ~~aut~~ ^{I.D.} autoclave, leave it in ^{I.D.} a water bath to cool down for 40 minutes.
 - a) ^{I.D.} Quick stir for 20 sec.
 - b) Pour plates by hand. Measure pH
 - c) Leave the plate on the bench to cool down.

(II) Buffer Peptone water 3 racks; 9 ml/tube.

- 1) Calculate the broth needed to prepare. $(9.5 \times 3 \times 75) + 3000 \text{ ml} = 25500 \text{ ml}$
Prepare 1.25 L x 2. could prepare all together in larger flask.
- 2) Calculate the ^{I.D.} broth needed. $(259 / 105 \text{ L})$ since ~~it~~ 20 g/L.
- 3) Marked the bottle where 1.25 L lies. Pour out DI water.
- 4) Weight 25 g using weighboats. carefully pour into the ~~flat~~ ^{I.D.} bottle.
- 5) Heat mix doesn't need heat sink ~~with~~
- 6) Use the dispenser to distribute 9 ml into each tube.
- 7) Autoclave 121°C 45 min lig.
- 8) Push the cap down, Leave the broth on the bench to cool down.

To Page No. _____

Witnessed & Understood by me,

yanon Deibany

Date ^{I.D.}

08/0

02/08/07

Invented by:

Uicaria Lopez

Recorded by:

Date

2/15/07

Project No. _____

Book No. _____

TITLE BHIA preparation.

From Page No. _____

Objective : To prepare 5 x 700 ml of Brain Heart Infusion Agar.

Procedure.

- 1) Marked the line where 700 ml lies. using cylinder and D₂O water.
- 2) Weight 96.4 g and pour into each flask.
- 3) Fill the ~~water~~ to 700 ml.
- 4) ~~Heat ^{pH}. Heat ^{pH} mix.~~ Cover with foil and heat mix.
- 5) Autoclave 121°C, 45 min, liq.
- 6) leave it in the waterbath (50°C), for 40 mins.
- 7) Pour plate. Measure pH.
- 8) leave it on the bench overnight to ~~cool down~~ ^{pH} solidify.

To Page No. _____

Witnessed & Understood by me,

yanam *[Signature]*

Date

02/13/07

Invented by:

[Signature]

Recorded by:

Date

2/15/07

Project No. _____

Book No. _____

TITLE BHIA preparation

From Page No. _____

Objective : To prepare 7 x 700 ml flasks of BHIA.

D.P. Procedure

- 1) Use cylinder to measure 700 ml DI water and use it to meet the volume of the flask at 700 ml.
- 2) Discard some of the water
- 3) Weight 26.4 g into each flask.
- 4) Fill in DI-water until it is 700 ml.
- 5) Drop a stir bar and cover with aluminium foil. Heat mix.
- 6) measure pH before autoclave.
- 7) Autoclave 121°C, 45 min, Lys.
- 8) Leave it in a waterbath (50°C) to cool down.
- 9) Pour plate. Measure pH of the media.
- 10) Leave it over night to solidify.

II) Weekly housekeeping.

- 1) Clean the waterbath.
- 2) Clean the balance and pH meter.

To Page No. _____

Witnessed & Understood by me,

Yaman Babur

Date

02/14/07

Invented by:

Recorded by:

Cristina Lopez

Date

2/15/07

Project No. _____

Book No. _____

TITLE Prepare 6 racks of 0.01% Peptone

From Page No. _____

Object: To prepare 6 racks of 0.01% Peptone. (9ml/tube)

1) Calculate the amount needed to prepare.

$$(9.5 \times 6 \times 75) + 300 = 4.6 \text{ L.}$$

2) Calculate the media (Peptone) needed.

$$1 \text{ L need } 0.1 \text{ g.}$$

$$4.6 \text{ L need } 0.46 \text{ g.}$$

3) ^{Y.D.} Use 1 L flask to prepare the media4) mark the glassware ^{Y.D.} at 4.6 L. by ^{Y.D.} using DI water and cylinder

5) Discard some of the water

6) Weight 0.46 g using small weight boat.

7) Fill the ^{Y.D.} DI water to the mark ^{Y.D.} 4.6 L.

8) Heat mix using magnetic bar to help mixing

9) Distribute to each tube by using dispenser

10) autoclave 121°C, 45 mins.

11) Leave them to cool down on the bench.

What else was done this day?

To Page No. _____

Witnessed & Understood by me,

Ganan Dubany

Date

02/15/07

Invented by:

Viviana Lynn

Recorded by:

Date

3/5/07

Project No. _____

Book No. _____

Prepare 7x700 ml of BHIA

LE _____

Page No. _____

Object: To prepare 7x700 of Brain Heart Infusion Agar.

- 1) ~~Measure 700 ml of DI~~ ^{Y.D}
- 2) Use Cylinder to measure 700 ml of DI water and use it to mark the glassware
- 3) ~~Weight 36.4 g of BHIA~~ ^{Y.D}
- 2) Discard some amount of water
- 3) Weight 36.4 g of BHIA and carefully pour the media into the flask.
- 4) Fill the water to the marked line.
- 5) Drop a magnetic bar and cover the flask with foil.
- 6) Heat mix til the agar ~~is~~ ^{Y.D} ~~and~~ ^{Y.D} homogenize with water. media homogenize.
- 7) Pour some of the media into a small container to measure the pH.
- 8) autoclave 121°C 45 min
- 9) Cool down in ~~the water~~ ^{Y.D} 45°C water bath. for 40 mins.
- 10) Stir for 20 seconds.
- 11) Pour plate
- 12) leave it to cool down on the bench.

Do daily housekeeping for lab 344 & Clean ~~lab~~ ^{Y.D} ~~lab~~
 Clean pH meter, balance, water bath & weekly housekeeping lab 315
 Clean table in lab 305

To Page No. _____

Accessed & Understood by me,

Arman Dahiya

Date

02/20/07

Invented by:

Arman Dahiya

Recorded by:

Date

3/5/07

Project No. _____

Book No. _____

TITLE _____

Prepare 3 racks of 0.01% peptone

From Page No. _____

Objective: To prepare 3 racks of 0.01% peptone.

Procedures

- 1) 3 racks, the volume needed is $(3 \times 75 \times 9.5) + 300 = 2.5 \text{ L}$.
- 2) Weigh 0.25 g of peptone into the 4 L flask that marked the 2.5 L volume and contain some DI water (around 1.8 L)
- 3) After carefully pour the media, fill water to the volume of 2.5 L
- 4) heat mix
- 5) distribute to tubes by using automatic dispenser
- 6) Autoclave 121°C , 45 min
- 7) Leave on the bench to cool down.

what else was done this day?

To Page No. _____

Witnessed & Understood by me,

yanan Sabotung

Date

2/20/07

Invented by:

Juwana Lopez

Recorded by:

Date

3/5/07

TLE Prepare 10 x 700 ml BHI

Project No. _____

Book No. _____

om Page No. _____

Objective: To prepare 10 x 700 ml BHI.

- 1) Mark the 1 L flask by using cylinder and H_2O to mark the volume at 700 ml.
- 2) Discard some amount of water
- 3) Weigh 36.4 g of BHI
- 4) Carefully pour the agar into the flask
- 5) Fill the water into the volume marked.
- 6) Drop in magnetic bar and cover with foil. label ~~with~~^{to} on an autoclave tape.
- 7) Heat mix
- 8) Pour some amount into a beaker to measure pH
- 9) Autoclave 45 min 121°C .
- 10) Put it in ~~in~~ⁱⁿ a water bath (45°C) for 40 mins.
- 11) Mix 20 s.
- 12) Pour plate

To Page No. _____

Witnessed & Understood by me.

2/21/07 *[Signature]*

Date

2/21/07

Invented by:

[Signature]

Recorded by:

Date

Project No. _____

Book No. _____

TITLE Prepare 8x 700 ml of BHIA

From Page No. _____


Objective: To prepare 8x 700 ml BHIA

- 1) Use 750 ml Cylinder to measure 700 ml of DI water
- 2) Mark each flask with level of 700 ml. Discard some water
- 3) Weight 36.4g of BHIA (use 52g)
- 4) Carefully pour the media into the flask
- 5) ~~heat mix~~ and ~~cover~~ cover with foil and heat mix
- 6) Autoclave 121°C, 45 min
- 7) Cool in water bath (45°C) for 40 min.
- 8) Mix 20 sec.
- 9) Pour plate.

- Do house keeping for each lab.

To Page No. _____


Witnessed & Understood by me,



Date

2/22/07

Invented by:



Recorded by:

Date

3/5/07

Prepare BHI A 7 x 700 ml
 Prepare buffer phosphate buffer (BPB), 4 y.p., 2 rocks.

Project No. _____

Book No. _____

m Page No. _____

Objective: To prepare 4^{yp.} racks of buffer phosphate buffer.

- 1) Prepare ~~(2 x 2.5 x 75) + 300~~^{yp.} = 425 L of BPB (4 x 9.5 x 75) + 300 ≈ 4 L.
- 2) Bring out the phosphate buffer stock solution
- 3) Since 1.25 ml ; 1 L need 5 ml
- 4) Use cylinder^L to measure ~~5~~ ^{yp.} ml DI-water and use it to mark 4 L of 6L flask
- 5) Use small cylinder to measure 5 ml ^{yp.} and fill in the flask
- 6) Heat mix
- 7) Distribute to each tube by using dispenser.
- 8) Autoclave 121°C, 45 min.
- 9) Cool on bench.

II Prepare BHI A 7 x 700 ml

- 1) Mark the point of 700 ml by using cylinder and DI water
- 2) Pour out some amount of water.
- 3) Weigh 86.4 g of BHI A and carefully pour into the flask
- 4) Fill the water to 700 ml
- 5) Drop magnetic bar and cover with flask foil^{yp.}.
- 6) Label on an autoclave tape a batch letter.
- 7) Heat mix
- 8) Autoclave 121°C, 45 min
- 9) Cool in water bath (48°C) 40 min.
- 10) Pour plate.

- Do housekeeping for lab 314, 314, 305

To Page No. _____

Witnessed & Understood by me,

Amrutha

Date

2/27/07

Invented by:

Mariana Kozu

Date

3/5/07

Recorded by:

Project No. _____

Book No. _____

TITLE Prepare 2x 700 TST, 2x 700 MOX, 4 racks of BPW.

From Page No. _____

- Objective:
- ① Prepare 2x 700 ml TST
 - ② Prepare 2x 700 ml MOX
 - ③ Prepare 4 racks 9.5 ml of Buffer peptone water.

Procedures.

I Prepare 2x 700 ml TST.

- 1) Use cylinder to weight DI water 700 ml and pour the water into the flask.
- 2) Mark the glass where the water lies.
- 3) Discard some water.
- 4) Use weight boat to weight 28 g of TSA media (1 L use 40
0.7 L use 28)
- 5) Carefully pour the media into the flask.
- 6) Pour DI water to the marked line.
- 7) Drop a stir bar and cover with foil
- 8) heat & thoroughly mix
- 9) Pour some media into ¹⁰ for pH measurement.
- 10) Sterilize by autoclave at 121°C, 45 min
- 11) Cool down in water bath (48°C) for 40 mins.
- 12) Pour plate

II Prepare 2x 700 ml MOX (oxford medium base)


- 1) Use cylinder to weight DI water 700 ml and pour the water into the flask
- 2) Mark the water line and discard some water
- 3) Weight 40.29 g of media on a weight boat (57.59 / 1 L , 0.7 L → 40.25)
- 4) Drop in a magnetic bar, ~~and~~ cover with foil, heat mix
- 5) Sterilize by autoclaving at 121°C, 45 min.
- 6) Cool down in water bath (48°C)
- 7) Prepare modified oxford medium supplement by rehydrate the supplement by adding 10 ml of DI water.
- 8) Add supplement to the ~~100~~ media after the media cool down (10 ml / L = 7 ml in
- 9) Pour plate
- 10) Leave to cool down on the bench.

IV. Prepare 4 racks of Buffered peptone water (9 ml)

- 1) Calculate the amount of media (9.5 × 4 × 75) + 300 = 3.2 L.
- 2) Marked 3.2 L ~~water~~ line by using cylinder and DI water
- 3) Weight 64 g of BPW (20 g / 1 L of water)

To Page No. _____


Witnessed & Understood by me,



Date

3/7/07

Invented by:



Recorded by:

Date


3/26/07

E _____

Page No. 20

- 4.) heat mix, ~~to~~ ^{to} ~~disp.~~
- 5.) Distribute into tubes by using dispenser.
- 6.) Autoclave 121°C, 45 min.
- 7.) leave on a bench to cool down.

Do daily housekeeping in lab ~~315~~ ^{41D} 315
 Do weekly housekeeping in lab 314, 315, 305
 Do inventory check for lab 314, 315 equipment.



To Page No. _____

Reviewed & Understood by me,
aman Pathanj

Date
3/7/07

Invented by:
Wendy Joyce
Recorded by:

Date
5/26/07

Project No. _____

Book No. _____

TITLE Prepare 9x700ml TSA, 8x1L PBS, 5x500ml 1x PBS

From Page No. _____

Objective: To prepare 9x700 ml of TSA. (3.1488 A-I)
 To prepare 8x 1L of 1x PBS and 5x 500ml of 1x PBS

I) Prepare 9x700 ml of TSA

- 1) Use cylinder to measure 700 ml of DI water and pour it into 1L flask.
- 2) Mark the line of water and pour out some water.
- 3) Weigh 28 g of media using weighboat and balance. (40 g/L)
- 4) Carefully pour the TSA into the flask. Add water to 700 ml.
- 5) Heat mix using stir bar, also cover with foil.
- 6) Autoclave 121°C, 45 min.
- 7) Cool down in water bath (48°C) for 40 min.
- 8) Per plate

II 1x PBS. (1L and 500 ml)

- 1) For 1L 1x PBS measure 900 ml of DI water and autoclave 121°C, 45 min
 For 500ml 1x PBS measure 450 ml of DI water and autoclave 121°C, 45 min.
- 2) Cool down the sterile water on the bench.
- 3) Add 100 ml of sterile 10x PBS into 900 ml sterile water to make 1x PBS 1 litre
 Add 50 ml of sterile 10x PBS into 450 ml sterile water to make 1x PBS 500ml
- 4) Cap tightly and shake vigorously to thoroughly mix it.

To Page No. _____

Witnessed & Understood by me

Gauri Bhatnagar

Date

3/7/07

Invented by:

Chiranjeev Singh

Recorded by:

Date

3/7/07

Project No. _____

Book No. _____

E Prepare 9 x 700 ml of XLD.

Page No. _____

Objective: Prepare 9 x 700 ml of XLD. (3/6/07)

Procedure

- 1) Use cylinder to measure 700 ml of DI water and pour it into 1L flask
- 2) ~~Discard~~ mark the line of the water and discard some amount of water
- 3) Weigh 38.5 g of XLD agar (55g/L) and carefully pour it into the flask
- 4) Pour DI water to ~~the~~ 700 ml
- 5) Drop stir bar and cover with foil, label
- 6) heat mix until media boils.
- 7) Cool down on the water bath (50°C) 40 mins.
- 8) Pour plate

Do housekeeping for the lab 315.
Continue with the inventory check of the equipment

To Page No. _____

Observed & Understood by me,

Arden

Date

3/12/07

Invented by:

Arden

Recorded by:

Date

3/24/01

Project No. _____

Book No. _____

TITLE Prepare 1x700 TSA, 1x700 Max, 2x9 ml x 2 racks BPN. (3)

From Page No. _____

Objective. To prepare 1x700 ml TSA, 1x700 ml Max, ^{2x} 9 ml x 2 racks of Buffer peptone water

Procedures.

I Prepare 1x700 ml of TSA

- 1) Use cylinder to measure 900 ml DI water and mark the water line
- 2) Discard some water
- 3) Weight 28 g of TSA and pour into the flask (1 L \rightarrow 40g)
- 4) Add water to the marked line 0.2 L \rightarrow 28 g
- 5) Heat mix using stir bar while covering with foil.
- 6) Autoclave, liq. 121°C 45 mins
- 7) Cool down in water bath (50°C) for 40 mins.
- 8) Pour plate using aseptic technique

II. Prepare 1x700 ml Max.

- 1) Mark the flask at the volume of 700 ml using DI water and cylinder. to measure the amount of water
- 2) Discard some water
- 3) Weight 40.25 g ^{1/2} of oxford medium base and pour it into a flask (0.25 g / 1 L, 0.2 L \rightarrow 40.25 g)
- 4) Add DI water to 700 ml
- 5) Autoclave, liq. 121°C 45 mins
- 6) Cool down in water bath for 40 mins (50°C)
- 7) Add 9 ml of modified oxford medium supplement
- 8) Mix 20 sec.
- 9) Pour plate

III Prepare 2 racks of BPN 9 ml.

- 1) Calculate the amount of media needed to prepare $(9.5 \times 25 \times 2) + 300 = 1.5 \text{ L}$
- 2) Use cylinder to measure 1.5 L DI water, marked the line.
- 3) Discard some water
- 4) Weight 30 g of BPN media (20 g/L) \therefore 1.5 L you need 30 g
- 5) Drop in magnetic bar and heat mix while covers the flask with foil
- 6) Cool down on bench
- 7) Distribute into tubes. using dispenser.
- 8) Autoclave liq. 121°C, 45 mins.
- 9) Cool down on bench.

DIP - Collect the hazard disposal.


To Page No. _____

Witnessed & Understood by me,

Date

Invented by:

Date



3/12/07

Recorded by:

3/24/07

Project No. _____

Book No. _____

Prepare 15 x 100 ml of TSA

Page No. _____

Objective: To prepare 15 ~~flask~~⁹¹⁰ 15 x 100 ml of TSA

Procedures

Use Cylinder to measure 100 ml of DI-water, Use marker to mark the line. Discard any water

Weight 22 g of TSA media (40g/L)

Add water to 100 ml

Heat mix

Autoclave 121°C, 45 mins.

Cool down in a water bath (50°C) 40 mins.

Pour plate

Do housekeeping for labs 314, 315, 305.

Collects kills

To Page No. _____

Read & Understood by me,

Navdeep Singh

Date

3/12/07

Invented by:

Navdeep Singh

Date

3/24/07

Recorded by:

Project No. _____

Book No. _____

TITLE Prepare 2x 100 ml of XLD

From Page No. _____

Objective: To prepare 2x 100 ml of XLD.

Procedure.

- 1) Use cylinder to measure 100 ml of DI water. Use label to mark the line
- 2) Discard some amount of water
- 3) Weigh 38.5 g of XLD media and carefully pour into the flask. (55g/L)
- 4) Turn on the water bath.
- 5) Heat mix ~~stir~~ until the media boils
- 6) Cool down in a water bath. 40 mins (50°C)
- 7) Pour plate

- Do the housekeeping for labs 314, 315, 305.
- Do kill run.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by:

Date

Janan Lami

3/13/07

Recorded by:

3/26/07

Project No. _____

Book No. _____

TITLE Prepare 2x700 ml TSA, Prepare 2x700 ml Max

From Page No. _____

Objective: To prepare 2x700 ml TSA and 2x700 ml Max

Procedure. I) 2x700 ml TSA

1. Use cylinder to measure 900 ml of DI water.
2. Mark the water line.
3. Pour out some water.
4. Weigh 28 g of TSA media (40 g/L) using weight boat.
5. Carefully pour the media into the flask, add DI water to 700 ml
6. Drop stir bar and cover with foil.
7. Heat mix
8. Auto clave, Liq 15 min.
9. Cool down in a water bath 40 mins (50°C)
10. Pour plate

II 2x700 Max

- 1) Use cylinder to measure 700 ml of DI water, use marker to mark the volume.
- 2) ~~10~~ Discard some water
- 3) Weigh 40.25 of media on a weight boat (57.5 g/L)
- 4) Add DI water to 700 ml
- 5) Auto clave, Liq 18°C 45 mins.
- 6) Cool down in water bath for 40 mins (50°C)
- 7) Add 1 ml of modified oxford medium supplement in to the ^{Y.D} agar flask.
- 8) Mix for 20 sec.
- 9) Pour plate

- Do the housekeeping for lab 314, 315, 305
- collect hazard kits

To Page No. _____

Witnessed & Understood by me,




Date

3/28/07

Invented by:

Recorded by:



Date

4/5/07

Project No. _____

Book No. _____

TITLE To Prepare 4 racks of BPH.

From Page No. _____

Objective: To prepare 9ml of buffer peptone water

1. Calculate the ~~media~~^{PO} amount of media that is needed to prepare $(9 \times 7.5 - 9.5 \times 75 \times 4) \div 300$
2. Use Cylinder to measure 3.2 L water and mark where the water lies.
3. Discard some amount of water.
4. Using weight boat to ~~measure~~^{weight} 64 g of BPH (20 g/L)
5. Carefully pour the media into the flask
6. Add DI water to 3.2 L.
7. heat mix until the media dissolve. - broths don't need to be heated
8. Distribute into tubes using dispenser.
9. Autoclave 121°C 45 min
10. Cool down on bench.

- Do housekeeping room 3/4

- Do kill run.

- Do inventory check of the expired media.

To Page No. _____

Witnessed & Understood by me,

Yamir Duly

Date

3/28/07

Invented by:

Urvash Jaisa

Recorded by:

Date

4/5/07

TITLE To prepare 15 x 900 ml TSA, Prepare ^{3 racks of} 10% sodium thiosulfate ~~3 racks of~~ Project No. _____
Book No. _____

From Page No. _____

Objective 2) To prepare 15 x 900 ml TSA

- 1) Use cylinder to measure 900 ml DI water, use marker to mark the volume
- 2) Discard some water
- 3) Weight 28 g. of TSA media (40 g / L)
- 4) Carefully pour the media in to the flask
- 5) Add DI water to 900 ml
- 6) Autoclave 121°C 45 min.
- 7) Cool in 50°C water bath. for 45 min.
- 8) Pour plate
- 9) Cool down on bench.

II Prepare 10% sodium thiosulfate 3 racks. (9 ml)

- 1) Use ⁹⁰⁰ cylinder calculate the amount of media needed. $(3 \times 9.5 \times 9.5) + 300 = 2.5 \text{ L.}$
- 2) Use cylinder to measure the volume and mark the flask.
- 3) Weight 250 grams of sodium thiosulfate (100 g / L.) and add it into flask.
- 4) beat 90 Mix thoroughly.
- 5) Distribute into tubes using dispenser.
- 6) Autoclave 121°C 45 min
- 7) Cool down on bench.

8) Do housekeeping in lab 315

To Page No. _____

Witnessed & Understood by me,

Jan Rij

Date

4/5/07

Invented by:

Recorded by:

Vivian Jaya

Date

4/5/07

Project No. _____

Book No. _____

TITLE ~~Report~~, 3 rocks of BPN.

From Page No. _____

- 9/28/07 :- The media is prepared but the pressure of the autoclave has not down to zero.
- Do the housekeeping for room 314, 315, 305

Should have informed me of the situation

To Page No. _____

Witnessed & Understood by me,

gan Duj

Date

9/5/07

Invented by:

Dwiana Leeza

Recorded by:

Date

9/5/07

TITLE Prepare 10xTSA (700 ml)

From Page No. _____

Objective: To prepare 10x700ml TSA

Procedures-

- 1) Use cylinder to measure 700 ml of D₂O water, pour into flask, and mark the volume.
- 2) Discard some water.
- 3) Weight 28 g of TSA into the flask (40g/L)
- 4) Add water to 700 ml
- 5) Heat mix while the flask is covered with foil.
- 6) Autoclave liq 45 min.
- 7) Cool down in ice water bath for 45 mins.
- 8) Pour plate.

- Do monthly housekeeping for lab 314
- Do daily housekeeping for lab 315 & 310

03-30-07 10:43 14

To Page No. _____

Witnessed & Understood by me,

Gyanou Rij

Date

4/5/07

Invented by:-

Sudana Singh

Recorded by:

Date

4/5/07

Project No. _____

Book No. _____

TITLE Prepare 10x700ml TSA.

From Page No. _____

Objective: To prepare 10x700ml TSA.

- 1) Mark the volume of 700 ml of the flask using DI water and cylinder.
- 2) Discard some amount of water.
- 3) Weigh 28 g of TSA and pour into the flask.
- 4) Add water to 700 ml.
- 5) Drop stir bar into the flask and cover with foil.
- 6) Heat mix.
- 7) Autoclave 121°C, 45 min.
- 8) Cool down in a 50°C water bath, 45 min.
- 9) Pour plate.
 - Do hazard disposal.
 - Do housekeeping for lab 314, 315, 305.

04-04-07 009:08 IN

TA

04-04-07 004:17 IN

Witnessed & Understood by me,

Yann Durr

Date

04/05/07

Invented by:

Recorded by:

Viviana Lopez

To Page No. _____

Date

4/5/07